

## MEDICAL EDUCATION

### Digital Breast Tomosynthesis versus Digital Mammography for Detection of Early-Stage Cancers Stratified by Grade: A TOSYMA Subanalysis

Weigel, S., Heindel, W., Decker, T., Weyer-Elberich, V., Kerschke, L., Gerf, J., Hence, H., for the TOSYMA Screening Trial Group  
Radiology 2023; 309(3):e231533. <https://doi.org/10.1148/radiol.231533>

#### Study Goal

- To compare invasive breast cancers detected with digital breast tomosynthesis (DBT) plus synthesized mammography (SM) versus digital mammography (DM) screening regarding tumor stage, histologic grade, patient age, and breast density.

#### Study Population/Method

- Multicenter, multivendor randomized controlled trial in the German mammography screening program
- 99,131 women imaged (49.9% DBT + SM, 50.1% DM screening)
- Median age: 57 years (range: 50 – 70 years)
- 83 readers of varying experience
- 55.2% of women presented with non-dense breasts, 44.8% presented with dense breasts
- >51% of images acquired with Hologic Mammography Systems
  - Including 11% high-resolution / Intelligent 2D images

#### Key Outcome Measures

- Compared invasive cancer detection rates (iCDRs), rate differences, and odds ratios between the study arms, stratified by stage, histological grade, age group, and breast density

#### Nottingham Scale

- Grade 2 and Grade 3 cancers are more aggressive cancers, and may have spread to lymph nodes or beyond

#### Results

The interim results of TOSYMA demonstrated that DBT+SM leads to:

- Increased detection of invasive cancers with DBT plus SM, as compared to DM alone
- Stage 1 iCDR difference: +21.6/1000 women 60-70 yrs w/dense breasts
- Stage 2+ iCDR difference: +2.3/1000 women 60-70 yrs w/dense breasts

#### Conclusion

- DBT plus SM appears to lead to higher detection of early-stage invasive breast cancers of Grade 2 or Grade 3 than DM screening, with the highest rate among women aged 60-70 years with dense breasts
- Higher detection rates of Grade 2 and Grade 3 invasive cancers with DBT plus SM, as compared to DM alone